

# Hipro Biotechnology Corp

Hipro is a high-tech enterprise based on international advanced medical technologies and brilliant self-innovations to provide first class products all over the world.

Hipro focuses on R&D, manufacturing, marketing, and relevant services of point-of-care products. Day in day out, Hipro is committed to innovate in the field of accurate diagnostic system including instrumentation, reagents, service, as well as data management. Distributed in 63 countries through a network of distributors, we can guarantee our customers a superior service and a product portfolio that meets the demand.

R&D Center: 3938 Trust Way, Hayward, CA94545, USA  
 Tel: +86 311 83832828  
 Service Hotline: 400-0191-606  
 Website: www.hipro.us

Email: [info@hipro.us](mailto:info@hipro.us)  
 Fax: +86 311 83859087



## A1

- Dimension: 410mm×280mm×355mm
- Weight: 5 Kg

## Test Channel

- Shake well automatically

## Large Display

- LCD Touch Screen
- Intelligent and intuitive operation interface

## Barcode Scanner

- Automatically reading of batch number and calibration information in the barcode

## SimplyRead

- One Step Operation
- Quantitative Sample collection
- Quantitative Reagent Filling
- Single Package For Liquid Reagents



## Printer

- Results printed automatically

## Build-in-Battery (optional)

- Charge 3Hr, run 36Hr non-stop

## Discharge port

- Tested reagents thrown out automatically

## 3-level Calibration System

- Hardware calibration
- Network calibration
- Third party calibration



### Collect sample



### Insert into cuvette



### Insert into channel



### Results printed



<b>SAS-Cov-2</b>	D-Dimer	HbA1c	PGII	<b>Tumor marker</b>
2019-nCoV IgM / IgG Antibody	HDL-C	Crea	HP	AFP
SAS-Cov-2-Antigen	TG	mALB	G17	CEA
<b>Infection</b>	BNP	Cys-C	<b>Liver function</b>	<b>Microelement</b>
CRP	NT-ProBNP	α1-MG	PA	25'-VD
hs-CRP	TCHO	β2-MG	TRF	Ferritin
PCT	FER	NGAL	ADA	BAP
SAA	HCY	RBP	<b>Autoimmunity</b>	GH
HBP	cTn I	<b>Rheumatism</b>	IgA	BGP
<b>Cardiovascular</b>	MYO	ASO	IgM	25-OH-VD
H-FABP	LDL-C	RF	IgE	
CK-MB	MPO	Anti-CCP	IgG	
LP (a)	cTnT	<b>Gastric Function</b>	<b>Hormone</b>	
Lp-PLA2	<b>Renal function</b>	PGI	HCG	